

CLOG ZAPPER INSTRUCTIONS:

1. Wash hands.
2. Open Clog Zapper package and take out supplies.
3. Open printed instruction paper to follow visual illustrations.
4. Turn feeding pump off and close roller clamp on administration set.
5. Insert clog zapper applicator into feeding tube side port (if present), if not, remove red adaptor from NG/NJ tube and insert directly into feeding tube.
6. Remove cap on 10cc syringe, insert end in cup of tap water, pull on plunger to fill syringe to 10cc mark on side of syringe, and put cap back on end.
7. **SHAKE! SHAKE! SHAKE!** until all powder is mixed and no lumps.
8. Remove cap from syringe.
9. Insert end of syringe into applicator, gently push end of plunger to instill 2-5cc of solution.
10. Remove applicator from tube and put plug in tube opening.
11. Keep applicator and syringe together.
12. Allow solution to stay in tubing 30-60 minutes, then
13. Fill 6cc syringe with tap water.
14. Remove plug from tube opening, insert end of 6cc syringe and push on plunger to flush tube.
15. If tube cannot be flushed, repeat steps #4 through #8.
16. After tube clears, dispose of supplies, and resume feeding.
17. If unable to get tube unclogged, call your doctor !

Note “CAUTION”, “STORAGE REQUIREMENTS”, and “INSTRUCTIONS FOR G-TUBE USE” at bottom of visual illustration sheet !!