



CVC LINE FLUSH

Flush unused CVC lines _____

Prepare work area.

Assemble supplies:

- Heparin filled syringe (one per each lumen)
- Alcohol swabs

Wash hands! Remove air from heparin syringe

Flush Line with Heparin :

1. Wipe end of valve with alcohol prep pad.
2. Connect Heparin syringe, and flush line, using a pulsing action.
3. Remove syringe.

CHANGE CLAVE VALVE WEEKLY:

1. Open valve package. Wipe end of valve with alcohol.
2. Attach prepared heparin syringe.
3. Prime valve with heparin.
4. Make sure catheter is clamped.
5. Wipe off junction of catheter and valve with alcohol prep pad..
6. Remove old valve. Screw on new valve.
7. Unclamp catheter.
8. Flush line as above. **DO NOT CLAMP CATHETER.**

Updated 05/06