

How to disconnect chemotherapy and begin hydration infusion:

- 1. When pump beeps “infusion complete” for chemotherapy, press the run/pause key then the on/off key to turn pump off.**
- 2. Disconnect chemo tubing from port, take out of pump and discard in chemotainer.**
- 3. Flush port with 10cc NS, checking for blood return.**
- 4. Spike bag upside down (since there is air in the bag) and squeeze bag to release air before fluid. When fluid reaches end of tubing, twist off oval yellow tab on flow stop to stop the flow.**
- 5. Load tubing in to pump and turn pump on by pushing on/off button.**
- 6. The pump will go through self checks after three beeps. The first screen will say “Alert! Perform routine maintenance” which will disappear quickly.**
- 7. When the pump arrives at the next screen it will have a black bar across the word “program”. Press yes.**
- 8. The third screen will have a black bar across “repeat rx”. Press yes.**
- 9. There will be two programs to choose from. Arrow down to “variable” for the hydration. Press yes.**
- 10. The pump will run through the program and stop at the “run” screen. Press “run” and it will start infusing.**