

ADULT HICKMAN FLUSH

Hickman lines should be flushed M-W-F using 3cc Heparin 100u/ml

Assemble supplies:

- Heparin-filled syringe (one per each lumen)
- Alcohol pads

Wash hands. Prepare work area.

Flush Line With Heparin:

- Remove air from syringe
- Wipe end of valve with alcohol prep pad
- Connect Heparin syringe, and flush line using pulsating motion
- Remove syringe

CHANGE MICROCLAVE VALVES WEEKLY

1. Make sure catheter is clamped
2. Wipe off junction of catheter and valve with alcohol prep pad
3. Have new valve ready
4. Remove old valve
5. Screw on new valve

mrh 10/04



Written 10/04
Revised
Reviewed