



IV INFUSION WITH INTERMATE AND SASH

Set up of Infusion

Assemble supplies:

- 2 Saline syringes Remove air
- 1 Heparin syringe Remove air
- Medicine in ***Intermate***
- Alcohol swabs
- Gloves

Wash hands. Prepare work area.

Flush Line with Saline

1. Scrub valve with alcohol.
2. Unscrew cap from saline syringe.
3. Screw saline syringe onto valve.
4. Inject Saline using push pause flush.
5. Unscrew saline flush.

Give Your Infusion

1. Screw end of ***Intermate*** tubing onto valve.
2. Open clamp on ***Intermate*** tubing and catheter tubing to begin infusion. Infusion should take approximately _____
3. When balloon is empty, clamp ***Intermate*** and remove Intermate from the valve.

Flush Line with Saline

1. Screw saline syringe onto valve.
2. Inject Saline using push pause method
2. Unscrew saline flush.

Flush Line with Heparin

1. Screw heparin syringe onto valve.
2. Inject heparin using push pause flush.
3. Unscrew heparin flush.

Flush Red line daily with Heparin.

