



PASV CVC LINE FLUSH

PASV CVC Flush

Flush unused CVC lines weekly

Prepare work area

Assemble supplies:

- Saline syringe (one for each lumen)
- Alcohol wipes
- Red cap

Wash hands!

Remove air from saline syringe

Flush line:

- 1. Open red cap and leave in plastic tray**
- 2. Wipe junction of cap and line with alcohol prep pad**
- 3. Remove cap and discard**
- 4. Connect saline syringe**
- 5. Flush using pulsing action**
- 6. Remove syringe and put on new red cap**